

Last updated: 23/01/2024 07:59:08

E. COLI/COLIFORM (25 TESTS)

Product group: **661** Product number: **777932**

The most basic test for bacterial contamination of a water supply is the test for total coliform bacteria. Total coliform counts give a general indication of the sanitary condition of a water supply.



For consistent results, an incubator should be used to culture the samples and the correct temperature.

(UV lamp supplied separately)

Product information

Colilert test is EPA approved method for total coliform and E.Coli in drinking water and E.Coli detection in wastewater.

Test Parameters: Coliform and E.Coli bacteria tests

Supplied with box containing pack of 25 Coli test bags and 25 ampules of growth media.

Used with:

- POTABLE WATER TEST KIT (778419)
- BALLAST CHECK 2 TEST KIT (777848)

Features

- Simple presence/absence test
- Easy to carry out on board ship
- Can be to test for Potable water and ballast water
- For Ballast water, dilute the sample by 1 part sample and 9 part deionised water (10 ml sample + 90 ml demin water = 100 ml sample)

Benefits

- Convenient method for measuring and monitoring bacteria and microbial activity in water,
- High quality and consistently reliable results

Specification

General

Invent Hazard Material (IMO/EU) classification	NA
--	----

Documents

[SDoC and MD for IHM](#)

Related products

Is accessory to

- 778419**
POTABLE WATER TEST KIT
- 777848**
BALLAST CHECK 2 TEST KIT

Is frequently bought together with

- 778483**
1 ML STERILE SYRINGES (30)
- 765018**
ENVIROCLEAN 25 LTR
- 777934**
CHLORINE DPD NO 1 (FREE) 250 TESTS
- 778414**
DILUTION POTS (40 PCS)
- 778413**
BACTERIA PLATES (40 PCS)

This page is printed from

<https://www.wilhelmsen.com/zh-hans/product-catalogue/products/marine-chemicals/test-kits-and-reagents/ballast-water-testing/e.coliform-test-set25-tests>